

MANUFACTURERS-DISTRIBUTORS OF PETROLEUM MINING & GENERAL LABORATORY TEST EQUIPMENT




(213)  
849-2179  
842-7179

CABLE: BICOBRAUN




BICO-BRAUN SHELL DESIGN ISCOMETER BATH A.S.T.M. D-445 NEW SOLID STATE MODULAR CONTROL PANEL




DULIN ROTAREX FOR ASPHALT TESTING TACHOMETER READ-OUT ELECTRIC BRAKE



UA BRAUN PULVERIZERS 1, 2, 3 H.P. V-BELT DRIVE




BRAUN CHIPMUNK CRUSHERS VD-400# HR. CAP. WD-800# HR. CAP.



BICO-BRAUN UD MODELS DIRECT DRIVE PULVERIZERS 2, 3 & 5 H.P.



MODEL # 13 PLUG WIRING DIAGRAM



BETWEEN RED & BLACK 240 VOLTS

BETWEEN RED & WHITE } 120 VOLTS  
AND BLACK & WHITE }