


BICO \* BRAUN INTERNATIONAL

MANUFACTURERS OF MINING, PETROLEUM AND ASSAY LAB EQUIPMENT SINCE 1888

BICO BRAUN  
ELECTRIC ASSAY FURNACE

STOCKING DISTRIBUTOR FOR DFC Ceramics  
FIRECLAY PRODUCTS AND UTENSILS

Catalog No. - 4151500EL 115v Furnace  
4242400EL 220v Furnace



\* FEATURES \*

- \* TEMPERATURE CONTROL
- \* DIRECT INDICATING PYROMETER
- \* 2300° F MAXIMUM TEMPERATURE
- \* CONVENIENT 115 volt 15 amp OR 220 volt 8 amp OPERATION
- \* INSIDE DIMENSIONS 9 X 10 X 6½"
- \* CAPACITY: FOUR TO SIX 30 Gram CRUCIBLES OR SIXTEEN TO TWENTY 1½" CUPELS
- \* APPROXIMATE SHIPPING WT. 180 lbs.
- \* SHIPPING DIMENSIONS 25½ X 26 X 22"

All Units Sold Complete Including Pyrometer, Temperature Control,  
Hearth Plate and Vent Hood.



BICO INCORPORATED  
PO BOX 6339  
BURBANK, CA 91510  
(213) 849-2179  
TELEX 194561 LSA Bico  
FAX # (818) 842-7976