





UA PULVERIZER



BICO
BRAUN INTERNATIONAL


ALLOW 1/2" AROUND BASE

Manufacturers of Petroleum and Mining Laboratory Test Equipment Since 1888





BICO, INC.
P.O. BOX 6339
BURBANK, CALIF. 91510, U.S.A. Cable BICOBRAUN

(213) 849-2179
(213) 842-7179



PULVERIZER SPEED 850 R.P.M.

ELEC. MOTOR SPEED
1725 R.P.M. 2HP FRAME
145 T - 3 PH. 220-140



BICO INC.
BURBANK CALIFORNIA
UA PULVERIZER

NOTE DO NOT TRY TO SCALE
SCHEMATIC DRAWING - NOT TO SCALE